

SoftForensics®



SCC Soft Computer's SoftForensics system – a flexible data solution for case management of Forensic Autopsy cases.

The SoftForensics system is a data and case management system to be used in Coroner and Forensic Autopsy and Testing facilities, by trained and experienced personnel. The SoftForensics system enables the user to manage and track the activities performed and the comprehensive data generated during the handling of forensic materials and reporting of cases associated with these materials, for the purpose of management of medico-legal autopsy cases.

The SoftForensics system supports the following functionalities:

- Documentation of decedent information and receipt in morgues, as well as subsequent autopsy procedures
- Tracking of legal documentation associated with the decedent and used for the determination of forensic vs. non-forensic case management
- Tracking of submitted evidence and release of such evidence to the appropriate entities
- Collation of test results, submitted evidence, and physical examinations into a comprehensive Forensic Case Report
- Audit trail of all activities performed for the case
- Security settings available to limit user access to functions in the system

FEATURES AND BENEFITS

Feature: Customizable case registration screens and information tabs

Benefit: Based on the case type – Homicide vs. Unknown vs. Anthropological vs. Non-Human Remains, the case registration screens may be defined to display only the needed data fields.

Feature: Dashboards configured to the logged in user role

Benefit: Based on the user role, the user will be provided a custom dashboard with only tiles relating to their job role.

Feature: Real time documentation of post mortem test requests

Benefit: The Pathologist is able to view the registration details for the case and available testing requisitions for real time ordering during the post mortem.

Feature: Confidentiality

Benefit: SoftForensics is an HTML based standalone application. It needs no integration with LIS or practice management systems as it contains its own case record database and setup files.

Feature: Inclusion of external reports and evidence in the case record

Benefit: SoftForensics provides the ability to store evidence and testing documents in the case record. Drag and drop boxes and scanning capability facilitate the addition of these documents.


Feature: Report generation capability

Benefit: SoftForensics provides the ability to generate autopsy reports using client custom templates and the MS Word Editor.

SAMPLE SCREENSHOTS

Custom Dashboard

Based on the logged in user role, the system displays the worklists and icons assigned to the role. Each icon denotes the number of cases qualifying for the particular worklist.



My PM's Worklist

Order additional testing real time during the post mortem. SoftForensics provides real time ordering to other SCC applications (if present) or print requisitions and labels for the specimens and tests ordered. The test catalogue and requisition template are user definable.

