

SoftWorkload[®]



SCC's SoftWorkload Measurement System and Reporting Frameworks are designed to report captured data in a standardized format. SoftWorkload units can be separated by department and procedure and can show work distribution between laboratories.

SoftWorkload is an SCC application that allows for the measurement and analysis of work performed across laboratory departments. SoftWorkload is designed to capture data from single site and multisite laboratories that use SCC applications. SoftWorkload is integrated with SoftLab[®], SoftMic[®], SoftBank[®], SoftCytogenetics[®], SoftPathDx[®], SoftTotalQC[®] and SoftMolecular[®]. Capturing of SoftWorkload data can be Ad-Hoc on demand and scheduled, automated from the crontab. SoftWorkload reports - weekly, monthly, quarterly, yearly - can be viewed, saved and/or printed, as well as scheduled to receive by email. Reporting formats include Adobe Reader PDF and Excel.


FEATURES AND BENEFITS

- | | |
|-----------------|---|
| Feature: | SoftWorkload allows using four methodologies, internationally recognized and utilized, to meet regulatory and non-regulatory needs. Regulatory standards can include: Canadian MIS Standards, American CAP- College of American Pathologists Standards, American LMIP - Laboratory Management Index Program Standards. A non-regulatory format can be used when the client's management is seeking to create their own internal reporting system. |
| Benefit: | Multiple reporting formats are available for use in one system. This gives users an enhanced and increased variety of reporting options to meet both regulatory and non-regulatory requirements. |
| <hr/> | |
| Feature: | SoftWorkload provides setup specifications to create Service Recipient (Patient) Activities and Non-Service Recipient (Non-Patient) Activities. |
| Benefit: | The all-inclusive capturing and reporting of Patient and Non-Patient Activities is accomplished all in one module. |
| <hr/> | |
| Feature: | With the Ad-Hoc Capturing feature, SoftWorkload users can run capturing for single or limited range of orders, units, Rx products, cases, equipment IDs, Lot records, and dates. Each module has its own dedicated Schedules of Capturing for capturing defined and saved activities. |
| Benefit: | Testing and validation of SoftWorkload setups is accomplished with ease utilizing SoftWorkload Ad-Hoc Capturing functionality in option CAPTURING. |
| <hr/> | |
| Feature: | SoftWorkload provides Setup Reports in PDF and Excel formats. |
| Benefit: | SoftWorkload Setup Reports (English and French, Excels and PDFs) can be viewed, saved and/or printed. They are very useful for testing and validation. |
| <hr/> | |
| Feature: | SoftWorkload capturing works from crontab on daily and monthly (extended) basis to capture data with the automated capturing feature. |
| Benefit: | Once implemented, SoftWorkload can be virtually hands-free, automatically capturing and reporting captured data. |
| <hr/> | |
| Feature: | WebScape SoftWorkload, modern look shows as a web page in a web browser. SoftWorkload Reports opens as separate tab in the same browser. |
| Benefit: | Clients can use SoftScape version or WebScape version, or both. |
| <hr/> | |
| Feature: | SoftWorkload integrated with SoftReports Print Scheduler to generate own on demand and scheduled reports in both PDF and Excel format. |
| Benefit: | SoftWorkload Statistical Reports can be scheduled to print or to receive by email. |
| <hr/> | |
| Feature: | Improved lookups for Tech ID, Workstation, Department, Location, etc. Narrow lookup search by single letter, number or character. |
| Benefit: | Shorter lookups, no lengthy scrolling to find needed data. |

SAMPLE SCREENSHOTS

SoftWorkload WebScape–Main Menu

An example of the Main Menu screen from new SCC SoftWorkload WebScape featuring tiles as options and sub options that system users can click on to navigate to defining screens.



SoftWorkload WebScape–SETUP– Schedule of Capturing

A screen example from new SCC SoftWorkload WebScape featuring Schedule of Capturing option from Setup.

SoftWorkload Excel Report–Selected Year, Total–Standard Report

An example of a yearly report run from option REPORTS, featuring PAT Patient Type as a COLUMN Category with SYS System and DEP Department as ROWS Sorting categories. SoftWorkload users have the option to run SoftReport PDFs and/or Excel Reports from SoftWorkload.

System	Department	Patient Type	E	I	O	TOTAL
L	SoftLab	MICN	0	4	0	4
L	SoftLab	VIR	0	4	0	4
M	SoftMic	QZD1	0	0	0	0
M	SoftMic	MICN	0	27	0	27
GRAND TOTAL:			0	41	0	41

SoftWorkload Statistical Reports – Parameters for Print Scheduler

A screen example showing parameters used to print SoftWorkload Statistical Reports from SoftReports Print Scheduler from SoftWorkload.