

# SoftMedia<sup>®</sup> for SoftBank<sup>®</sup>



SoftMedia is user-friendly and provides the tools necessary to store documents and digital images from a variety of sources in a central repository which can be easily linked to patient records.

SoftMedia is fully integrated with SoftBank to allow users to electronically store documents in SoftBank such as Requisitions or ABID forms as well as any other type of reference material or documents of interest.

As an integrated module, SoftMedia does not require any additional interfaces and runs on the SoftBank server itself. SoftMedia saves money by reducing or eliminating the need to store and retrieve paper records, which is costly in terms of both space and labor, all while promoting a paperless environment.

## FEATURES AND BENEFITS

**Feature:** Electronical storing of documents

**Benefit:** SoftMedia allows users to scan documents or add files to SoftBank orders. The system automatically attaches documents to the specific SoftBank order, without requiring the user to manually provide any information.

**Feature:** Easily viewable documents

**Benefit:** SoftMedia provides easy retrieval and preview of the document from the order for the user to review as needed.

**Feature:** Supporting multiple formats

**Benefit:** Patient documents and images can be stored from a wide variety of data sources including scanners, digital cameras and other imaging devices, instrument interfaces, e-mail, etc.

**Feature:** Batch scan large volumes of SoftBank documents

**Benefit:** SoftMedia's batch scan process allows users to quickly and efficiently scan large volumes of documents and automatically link the documents to the correct order number or case number. The scanned documents can be easily retrieved from other SCC applications for viewing purposes.

**Feature:** Import image options

**Benefit:** SoftMedia allows to store documents from instruments, utilizing the Import Utility to automatically attach to the correct order.

**Feature:** Workspaces

**Benefit:** SoftMedia provides separate libraries for documents, which allows a physical separation of documents, and can also be used to separate historical documents. This efficient and systematic feature improves work performance and streamlines the data storage process by providing maximum control over documents and multiple archive queries for specific documents and images. Some workspaces are compatible with external document management systems.


**Feature:** Manage document storage

**Benefit:** SoftMedia allows facilities to manage the storage life of SoftMedia stored documents, with automatic archiving and purging as specified by the facility.

## SAMPLE SCREENSHOTS


### Electronically stored Antibody Panel

Antibody Panel can be electronically stored and attached to a patient's antibody identification test order.



### Electronically stored Reference Lab Report

Reference Lab Reports can be stored and attached to a Patient Order, either through scanning the document or uploading the electronic file.



### SoftMedia Workspaces

SoftMedia workspaces can be created and used for different purposes. Documents such as Reagent Package Inserts could be organized and stored in a custom Shared Workspace and made easily accessible to all necessary users. Different views can be created to present the same workspaces differently to different users or for various purposes. Different perspectives can filter out files which don't need to be presented for certain functionality. Security can be added to Workspaces to allow access restriction.

